

MoleMax **HD** Family



Mole Max **HD**

Two touch screen monitors

Additional camera preview LCD screen

Maximum resolution by digital image transfer via HDMI

Total Body Mapping capability


Mole Max **HD** AMTBM

Advanced LED Flash with polarization filter

Fast and accurate automatic camera positioning

Auto mole count Software for fast lesion enumeration on TBM images

Patient positioning measurement



MoleMax **HD** modular

Flexible configuration and positioning of components at your needs

MoleMax HD camera and computer

Additional camera preview LCD screen

Portable Dermoscopy **HD** Imaging



DermDoc HR

Polarized ELM Illumination and HD Quality Imaging

Autofocus

Automatic mode selection, Micro&Macro

Proved MoleMax Software

SkinDOC

Fixed Focus full HD camera

Digital real time video

Easy self-installation

Proved MoleMax Software



PhotoMAX PRO DLFII

Polarized and non polarized mode

Live image display

Camera control by software

Wifi option

