


PhotoMAX PRO DLFII

NEW ADVANCED DIGITAL IMAGE CAPTURE



VERSATILE ADVANCED DIGITAL DERMOSCOPY SET

From clinical overview up to advanced Dermoscopy, all images in high resolution quality. Together with all new and redeveloped PhotoMax Documentation-Software, images are also displayed online on the Monitor.



System Requirements:

From Windows 7 and later, Graphic display adapter capable of 1024 x 768 or more with True Color (24 or 32 bit) resolution, one USB 2.0 port, 1 GB free hard drive space, Mouse, Keyboard

PhotoMAX PRO DLFII

Consists of:

- High resolution Canon SLR digital Camera
- New versatile ELM Lens including integrated rechargeable battery pack
- Canon Macro Lens for Clinical Imaging
- Unique new Camera software control
- Advanced documentation Software
- Charger
- Hard case transport box

Applications:

- Online Live image via USB
- Macro and clinical imaging capability
- ELM Interface for high resolution Dermoscopy
- Switchable between polarized and non-polarized mode
- Live image display on monitor & on the Camera
- Optical magnification x12 (on Screen x 30)
- Field view 9-20mm
- Diagnosis support using Expertizer
- Real time follow up examinations
- Networkable with any other MoleMax system or Software